

## Product Datasheet

### Bovine ANGPT1 ELISA Kit (orb339676)

|                          |                                             |
|--------------------------|---------------------------------------------|
| <b>Description</b>       | Bovine Angiopoietin-1 (ANGPT1)<br>ELISA kit |
| <b>Reactivity</b>        | Bovine                                      |
| <b>Range</b>             | 12.5 ng/mL-800 ng/mL                        |
| <b>Target</b>            | ANGPT1                                      |
| <b>Note</b>              | For research use only                       |
| <b>Application notes</b> | Detection Wavelength: 450 nm                |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates        |
| <b>Uniprot ID</b>        | <b>O18920</b>                               |
| <b>Sensitivity</b>       | 1.384 ng/mL                                 |
| <b>Expiration Date</b>   | Please enquire.                             |