

## Product Datasheet

### Chicken APOB ELISA Kit (orb339681)

|                          |                                              |
|--------------------------|----------------------------------------------|
| <b>Description</b>       | Chicken Apolipoprotein B (APOB)<br>ELISA kit |
| <b>Reactivity</b>        | Gallus                                       |
| <b>Range</b>             | 0.156 µg/mL-10 µg/mL                         |
| <b>Target</b>            | APOB                                         |
| <b>Note</b>              | For research use only                        |
| <b>Application notes</b> | Detection Wavelength: 450 nm                 |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates         |
| <b>Uniprot ID</b>        | <b>P11682</b>                                |
| <b>Sensitivity</b>       | 0.039 µg/mL                                  |
| <b>Expiration Date</b>   | Please enquire.                              |