

## Product Datasheet

### Human ADP ELISA Kit (orb339700)

|                          |                                                                      |
|--------------------------|----------------------------------------------------------------------|
| <b>Description</b>       | Human adiponectin, ADP ELISA Kit                                     |
| <b>Reactivity</b>        | Human                                                                |
| <b>Range</b>             | 1.562 ng/mL-100 ng/mL                                                |
| <b>Target</b>            | ADP                                                                  |
| <b>Note</b>              | For research use only                                                |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                         |
| <b>Sample Types</b>      | serum, plasma, cell culture supernates,<br>tissue homogenates, urine |
| <b>Uniprot ID</b>        | <b>Q15848</b>                                                        |
| <b>Sensitivity</b>       | 1.102 ng/mL                                                          |
| <b>Expiration Date</b>   | Please enquire.                                                      |