

Product Datasheet

Human Apo-A1 ELISA Kit (orb339707)

Description	Human apolipoprotein A1 (Apo-A1) ELISA Kit
Reactivity	Human
Range	1.4 ng/mL-1000 ng/mL
Target	APOA1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, saliva, urine
Uniprot ID	P02647
Sensitivity	0.7 ng/mL
Expiration Date	Please enquire.