

Product Datasheet

Human sTWEAK ELISA Kit (orb339721)

Description	Human soluble tumor necrosis factor-like weak inducer of apoptosis (sTWEAK) ELISA Kit
Reactivity	Human
Range	0.312 ng/ml-20 ng/ml
Target	TNFSF12
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O43508
Sensitivity	0.078 ng/ml
Expiration Date	Please enquire.