

Product Datasheet

Human TSP-1 ELISA Kit (orb339722)

Description	Human thrombospondin 1, TSP-1 ELISA Kit
Reactivity	Human
Range	2.7 ng/ml - 2000 ng/ml
Target	THBS1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, breast milk, saliva, tissue homogenates, cell culture supernates
Uniprot ID	P07996
Sensitivity	2.7 ng/ml
Expiration Date	Please enquire.