

Product Datasheet

Human TIMP-1 ELISA Kit (orb339725)

Description	Human Tissue Inhibitor Of Matrix Metalloprotease-1 (TIMP-1) ELISA KIT
Reactivity	Human
Range	0.39 ng/mL-25 ng/mL
Target	TIMP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P01033
Sensitivity	0.434 ng/mL
Expiration Date	Please enquire.