

Product Datasheet

Human VEGFR-2/Flk-1 ELISA Kit (orb339730)

Description	Human Vascular endothelial cell growth factor receptor 2, VEGFR-2/Flk-1 ELISA kit
Reactivity	Human
Range	0.39 ng/mL-25 ng/mL
Target	VEGFR-2/Flk-1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P35968
Sensitivity	0.098 ng/mL
Expiration Date	Please enquire.