

Product Datasheet

Human Vaspin ELISA Kit (orb339731)

Description	Human Visceral adipose-specific serine protease inhibitor, vaspin ELISA Kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	SERPINA12
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8IW75
Sensitivity	7.8 pg/mL
Expiration Date	Please enquire.