

Product Datasheet

Mouse ADP ELISA Kit (orb339732)

Description	Mouse adiponectin, ADP ELISA Kit
Reactivity	Mouse
Range	0.156 ng/ml-10 ng/ml.
Target	ADP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates.
Uniprot ID	Q60994
Sensitivity	0.039 ng/ml
Expiration Date	Please enquire.