

Product Datasheet

Mouse ANGPTL3 ELISA Kit (orb339734)

Description	Mouse Angiopoietin-related protein 3 (ANGPTL3) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Target	ANGPTL3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9R182
Sensitivity	10.24 pg/mL
Expiration Date	Please enquire.