

Product Datasheet

Mouse Apo-E ELISA Kit (orb339736)

Description	Mouse apolipoprotein E (Apo-E) ELISA Kit
Reactivity	Mouse
Range	6.25 ng/mL-400 ng/mL
Target	APOE
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P08226
Sensitivity	1.56 ng/mL
Expiration Date	Please enquire.