

Product Datasheet

Mouse TIMP-1 ELISA Kit (orb339746)

Description	Mouse tissue inhibitors of metalloproteinase 1, TIMP-1 ELISA Kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Target	TIMP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates
Uniprot ID	P12032
Sensitivity	0.110 ng/mL
Expiration Date	Please enquire.