

Product Datasheet

Rabbit Apo-B100 ELISA Kit (orb339750)

Description	Rabbit apolipoprotein B100 (Apo-B100) ELISA Kit
Reactivity	Rabbit
Range	78 ng/mL-5000 ng/mL
Target	Apo-B100
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P17165
Sensitivity	19.5 ng/mL
Expiration Date	Please enquire.