

Product Datasheet

Rat ICAM-1 ELISA Kit (orb339757)

Description	Rat intercellular adhesion molecule 1, ICAM-1 ELISA Kit
Reactivity	Rat
Range	78.125 pg/mL-5000 pg/mL
Target	ICAM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q00238
Sensitivity	24.321 pg/mL
Expiration Date	Please enquire.