

Product Datasheet

Human CEACAM1 ELISA Kit (orb339784)

Description	Human Carcinoembryonic antigen-related cell adhesion molecule 1 (CEACAM1) ELISA kit
Reactivity	Human
Range	0.78 ng/mL-50 ng/mL
Target	CEACAM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P13688
Sensitivity	0.195 ng/mL
Expiration Date	Please enquire.