

Product Datasheet

Human MCP1/MCAF ELISA Kit (orb339811)

Description	Human MCP1/MCAF ELISA kit
Reactivity	Human
Range	6.25 pg/ml-400 pg/ml.
Target	CCL2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, urine, saliva, cerebrospinal fluid (CSF), ascitic fluid, cell culture supernates, tissue homogenates.
Uniprot ID	P13500
Sensitivity	1.56 pg/ml.
Expiration Date	Please enquire.