

Product Datasheet

Human MPIF-1 ELISA Kit (orb339812)

Description	Human Myeloid Progenitor Inhibitory Factor 1, MPIF-1 ELISA Kit
Reactivity	Human
Range	39 pg/mL-2500 pg/mL
Target	CCL23
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P55773
Sensitivity	9.7 pg/mL
Expiration Date	Please enquire.