

Product Datasheet

Human TNFRSF17 ELISA Kit (orb339819)

Description	Human Tumor necrosis factor receptor superfamily member 17 (TNFRSF17) ELISA kit
Reactivity	Human
Range	7.8 pg/ml- 500 pg/ml
Target	TNFRSF17
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q02223
Sensitivity	1.95 pg/ml
Expiration Date	Please enquire.