

Product Datasheet

Mouse MCP-5 ELISA Kit (orb339828)

Description	Mouse Monocyte Chemotactic Protein 5 (MCP-5) ELISA Kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Target	MCP-5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q62401
Sensitivity	16.82 pg/mL
Expiration Date	Please enquire.