

Product Datasheet

Hamster adiponectin ELISA Kit (orb339856)

Description	Hamster adiponectin, C1Q and collagen domain containing (ADIPOQ) ELISA kit
Reactivity	Golden Hamster
Range	0.469 ng/mL-30 ng/mL
Target	ADIPOQ
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	C9DFG2
Sensitivity	0.117 ng/mL
Expiration Date	Please enquire.