

Product Datasheet

Human Lp-PLA2 ELISA Kit (orb339865)

Description	Human lipoprotein-associated phospholipase A2, Lp-PLA2 ELISA Kit
Reactivity	Human
Range	1.56 ng/mL-100 ng/mL
Target	PLA2G7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q13093
Sensitivity	0.39 ng/mL
Expiration Date	Please enquire.