

Product Datasheet

Human NGAL/MMP-9 complex ELISA Kit (orb339869)

Description	Human NGAL/MMP-9 complex ELISA kit
Reactivity	Human
Range	0.313 ng/mL-20 ng/mL
Target	NGAL/MMP-9 complex
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, saliva, cell culture supernates, tissue homogenates
Sensitivity	0.029 ng/mL
Expiration Date	Please enquire.