

## Product Datasheet

### Horse Rotavirus ELISA Kit (orb339911)

|                          |                                                |
|--------------------------|------------------------------------------------|
| <b>Description</b>       | Horse rotavirus (RV) antigen<br>(Ag) ELISA kit |
| <b>Reactivity</b>        | Equine                                         |
| <b>Range</b>             | Request Information                            |
| <b>Target</b>            | RV Ag                                          |
| <b>Note</b>              | For research use only                          |
| <b>Application notes</b> | Detection Wavelength: 450 nm                   |
| <b>Sample Types</b>      | feces                                          |
| <b>Sensitivity</b>       | Request Information                            |
| <b>Expiration Date</b>   | Please enquire.                                |