

Product Datasheet

Human CSTB/CST6/STFB ELISA Kit (orb340024)

Description	Human Cystatin-B (CSTB/CST6/STFB) ELISA kit
Reactivity	Human
Range	0.313 ng/mL-20 ng/mL
Target	CSTB
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, urine, saliva, tissue homogenates, cell lysates
Uniprot ID	P04080
Sensitivity	0.233 ng/mL
Expiration Date	Please enquire.