



Product Datasheet

NPTX2 antibody (orb34003)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	NPTX2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47042
Uniprot ID	P47972
NCBI	NP_002514.1
Dilution Range	WB 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human te...