

Product Datasheet

Human TNFRSF9 ELISA Kit (orb340032)

Description	Human Tumor necrosis factor receptor superfamily member 9 (TNFRSF9) ELISA kit
Reactivity	Human
Range	7.8 pg/mL-500 pg/mL
Target	TNFRSF9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q07011
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.