



Product Datasheet

Olfm1 antibody (orb34004)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	OLFM1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55343
Uniprot ID	Q99784
NCBI	NP_055094.1 , NP_006325.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of Olfm1 antibody