

## Product Datasheet

### Fish Testosterone ELISA Kit (orb340065)

|                          |                                   |
|--------------------------|-----------------------------------|
| <b>Description</b>       | Fish Testosterone (T) ELISA Kit   |
| <b>Reactivity</b>        | Fish                              |
| <b>Range</b>             | 0.1 ng/mL-20 ng/mL                |
| <b>Target</b>            | T                                 |
| <b>Note</b>              | For research use only             |
| <b>Application notes</b> | Detection Wavelength: 450 nm      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates |
| <b>Sensitivity</b>       | 0.05 ng/mL                        |
| <b>Expiration Date</b>   | Please enquire.                   |