

## Product Datasheet

### Human Chitinase-3-like Protein 1 ELISA Kit (orb340082)

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Description</b>       | Human Chitinase-3-like Protein 1<br>(YKL-40/CHI3L1) ELISA Kit |
| <b>Reactivity</b>        | Human                                                         |
| <b>Range</b>             | 46.875 pg/mL-3000 pg/mL                                       |
| <b>Target</b>            | CHI3L1                                                        |
| <b>Note</b>              | For research use only                                         |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                  |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates                          |
| <b>Uniprot ID</b>        | <b>P36222</b>                                                 |
| <b>Sensitivity</b>       | 78.065 pg/mL                                                  |
| <b>Expiration Date</b>   | Please enquire.                                               |