

## Product Datasheet

### Monkey Progesterone ELISA Kit (orb340103)

|                          |                                      |
|--------------------------|--------------------------------------|
| <b>Description</b>       | Monkey Progesterone (PROG) ELISA Kit |
| <b>Reactivity</b>        | Monkey                               |
| <b>Range</b>             | 0.15 ng/ml-70 ng/ml.                 |
| <b>Target</b>            | PROG                                 |
| <b>Note</b>              | For research use only                |
| <b>Application notes</b> | Detection Wavelength: 450 nm         |
| <b>Sample Types</b>      | serum, plasma                        |
| <b>Sensitivity</b>       | 0.2 ng/mL                            |
| <b>Expiration Date</b>   | Please enquire.                      |