

Product Datasheet

CD133 antibody (orb340142)

Description	CD133 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD133. Please contact us for the exact immunogen sequence. The peptide is available as orb374956.
Target	CD133. CD133 antibody detects endogenous levels of total CD133
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	O54990, O43490
NCBI	001139319, 001145847, 11
Entrez	8842
Dilution Range	IHC-P:1:200, WB: 1:500, IF/ICC: 1:200
Expiration Date	12 months from date of receipt.