

Product Datasheet

ZO-1 antibody (orb340146)

Description

ZO-1 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P

Applications

Immunogen

KLH conjugated synthetic peptide derived from human ZO-1. Please contact us for the exact immunogen sequence. The peptide is available as orb12904.

Target

ZO-1. Reconnized detects endogenous levels of total ZO-1 antibody

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

Uniprot ID

[P39447](#), [Q07157](#)

NCBI

[003248](#), [003257](#)

Entrez


[7082](#)

Dilution Range


IHC-P:1:100-500

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of mouse ut...



IHC-P image of rat colon tissue using an...