




Product Datasheet

Vimentin antibody (orb34038)

Description	Rabbit polyclonal antibody to Vimentin
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	VIM
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53652
Uniprot ID	P08670
NCBI	NP_003371.2
Dilution Range	IF/ICC: 1:200, WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell line ...



Immunofluorescence analysis of SY5Y cell...