

Product Datasheet

Caspase 3 antibody (orb340700)

Description

Caspase 3 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Caspase 3. Please contact us for the exact immunogen sequence. The peptide is available as orb374962.

Target

Caspase 3. Caspase 3 antibody detects endogenous levels of total Caspase 3

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

32 kDa

Uniprot ID

P70677, **P42574**, **P55213**

Entrez


836

Dilution Range


IHC-P: 1:100-500, WB: 1:100-2000

Expiration Date


12 months from date of receipt.



Western blot analysis of rat ovary (lane...)



Immunohistochemical staining of paraffin...



Immunohistochemical staining of rat blad...