




Product Datasheet

Collagen IX antibody (orb34075)

Description	Rabbit polyclonal antibody to Collagen IX
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	COL9A1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	91869
Uniprot ID	P20849
NCBI	NP_001842.3 , NP_511040.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 cell line...



Western blot analysis of MCF-7 cell line...