



Product Datasheet

OGT antibody (orb340807)

Description	Rabbit polyclonal antibody to OGT
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, IP, WB
Immunogen	Recombinant full length protein of human OGT
Target	OGT
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	O15294 , P56558 , Q8CGY8
Entrez	26295 , 8473 , 108155
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:50-1:200, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of OGT expression ...



Immunohistochemical analysis of OGT stai...



Immunofluorescent analysis of OGT staini...