



Product Datasheet

CPT2 antibody (orb340854)


Description	Rabbit polyclonal antibody to CPT2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide of human CPT2
Target	CPT2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P23786 , P18886 , P52825
Entrez	25413 , 12896 , 1376
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IF/ICC: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of CPT2 expression...



Immunohistochemical analysis of CPT2 sta...



Immunofluorescent analysis of CPT2 stain...