



Product Datasheet

PRDM14 antibody (orb340932)


Description	Rabbit polyclonal antibody to PRDM14
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IH, IP, WB
Immunogen	Recombinant full length protein of human PRDM14
Target	PRDM14
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9GZV8 , E9Q3T6
Entrez	63978 , 383491
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:50-1:200, IP: 1:20-1:100, ChIP: 1:20-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of PRDM14 expressi...



Immunohistochemical analysis of PRDM14 S...



Immunofluorescent analysis of PRDM14 sta...