

Product Datasheet

LDHD antibody (Biotin) (orb34094)

Description	LDHD antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Probable D-lactate dehydrogenase, mitochondrial protein (310-507AA)
Target	LDHD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86WU2
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.