



Product Datasheet

ELAVL3 antibody (orb340994)


Description	Rabbit polyclonal antibody to ELAVL3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, IP, WB
Immunogen	KLH-conjugated synthetic peptide of human HuC
Target	ELAVL3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q60900, Q14576
Entrez	15571, 1995
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HuC expression ...



Immunohistochemical analysis of HuC stai...



Immunofluorescent analysis of HuC staini...