


Product Datasheet

RAD54L2 antibody (orb341001)

Description	Rabbit polyclonal antibody to RAD54L2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant full length protein of human RAD54L2
Target	RAD54L2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9Y4B4
Entrez	23132
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:20-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of RAD54L2 express...