



Product Datasheet

COG2 antibody (orb341018)

Description	Rabbit polyclonal antibody to COG2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	Recombinant full length protein of human COG2
Target	COG2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q921L5 , Q14746
Entrez	22796 , 76332
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:20-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of COG2 expression...



Immunohistochemical analysis of COG2 sta...