



Product Datasheet

CACNA2D4 antibody (orb34110)


Description	CACNA2D4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Voltage-dependent calcium channel subunit alpha-2/delta-4 protein (934-1115AA)
Target	CACNA2D4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z3S7
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human sk...



Immunofluorescence analysis of HepG2 cel...



Immunohistochemical staining of human he...