


Product Datasheet

GTF2H3 antibody (orb341249)

Description	Rabbit polyclonal antibody to GTF2H3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant full length protein of human TFIIH p34
Target	GTF2H3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8VD76 , Q13889
Entrez	2967 , 209357
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of TFIIH p34 expression.