


Product Datasheet

FTL antibody (Biotin) (orb341492)

Description	Goat polyclonal antibody against FTL conjugated to Biotin
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	FTL, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 20kDa band observed in Rat Brain lysates (calculated MW of 20.7kDa according to Rat NP_071945.3). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 1-3µg/ml.
Protein Sequence	GEYLFERLTLKHD.
Clonality	Polyclonal
MW	20
Entrez	2512
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of Rat...