

Product Datasheet

SCARB2 antibody (Biotin) (orb341498)

Description	Goat polyclonal antibody against SCARB2 conjugated to Biotin
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	LIMP2 / SCARB2, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 80kDa band observed in Human Brain (Frontal Cortex) lysates (calculated MW of 54.3kDa according to NP_005497.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Fujita et al, Biochem Biophys Re
Protein Sequence	NKANIQFGDNGTTIS.
Clonality	Polyclonal
MW	54.3
Entrez	950
Dilution Range	ELISA: 1:4000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of Hum...