


Product Datasheet

Adig antibody (orb341504)

Description	Goat polyclonal antibody to Adig
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	Adipogenin (mouse, aa41-52)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments in Mouse Adipose lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).. IHC: Paraffin embedded Human Testis. Recommended concentration: 5µg/ml.
Protein Sequence	SQESKEDDSDLC
Clonality	Polyclonal
MW	9.36
Dilution Range	ELISA: 1:2000, WB: 1 µg/ml), IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Testis tissue showing Adig antibody binding.