


Product Datasheet

NIBR-17 (orb341526)

Description	NIBR-17 is a pan class I PI3K inhibitor with IC50 values of 1, 9.2, 9 and 20 nM respectively for PI3...
Note	For research use only
Formula	C ₁₈ H ₂₀ N ₈ O ₂
MW	380.4
CAS Number	[944396-88-7]
SMILES	<chem>COC1=NC=C(C=C1)NC2=NC(=NC(=C2)C3=CN=C(N=C3)N)N4CCOCC4</chem>
Expiration Date	12 months from date of receipt.



Chemical
structure of
NIBR-17