

Product Datasheet

TRAPPC6B antibody (FITC) (orb34153)

Description	TRAPPC6B antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Trafficking protein particle complex subunit 6B protein (1-130AA)
Target	TRAPPC6B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86SZ2
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.